






BACK-STATION-COMBI






Isomed2000 Back-Station-Combi

- kombinira dva sustava u jednom uređaju s vlastitim pogonom dinamometra
- vlastiti sustav upravljanja
- vlastiti nadzor nad glavom dinamometra
- brzo pomicanje i klikanje priključka dinamometra

Left / Right side back rotation



Back flexion extension



Dinamometar povezan s jedinicom za rotaciju trupa



Dinamometar povezan s jedinicom za fleksiju/ekstenziju trupa



Pomaknite i kliknite glavu dinamometra



Dinamometar spojen sa jedinicom za rotaciju trupa



BACK-STATION-COMBI

TEHNIČKI PODACI

FLEKSIJA/EKSTENZIJA TRUPA

Back Flexion/Extension		
Max obrtni moment concentric	1000Nm	
Max obrtni moment eccentric	1000Nm	
Max brzina	200°/s	

Back Flexion/Extension electrical adjustments		
Seat-Height	up	down
Seat-Displacement	forward	backwards
Food-plate height	up	down
Pelvis pad displacement	forward	backwards
Pelvis pad height	up	down



BACK-STATION-COMBI

TEHNIČKI PODACI

ROTACIJA TRUPA

Max obrtni moment concentric	750Nm	
Max obrtni moment eccentric	750Nm	
Max brzina	180°/s	

Electrical adjustments		
Shine-bone fixation	forward	backward
food-plate height	up	down

Mechanical adjustments		
Side-thigh pad wide	narrow	wide
thigh pad angle		
Hip pad displacement	forward	backward
Shoulder front pad height	up	down
Shoulder front pad fixation	upward	downward

BACK-STATION-COMBI

2 D DRAWING

