








## Technical data


	Gear 1	Gear 2
Speed, maximum	0,8 m/s	1,2 m/s
Force with 750/560 Dynamometer (IsoMed E750-560)	8850 N	5850 N
Force with 700/560 Dynamometer (IsoMed E700-560)	8250 N	5500 N
Acceleration, max. adjustable in 0,02 m/s <sup>2</sup> steps*	7,5 m/s <sup>2</sup>	7,5 m/s²
Deceleration in 0,02 m/s <sup>2</sup> steps	10,0 m/s²	10,0 m/s²
Displacement of sled, maximum	0-51 cm	0-76 cm
Weight of Legpress athletic	485 KG	



## Dimensions







All measurements include the EASY MOVE wheels. Please be aware that the IsoMed 2000 may need additional space and you may need additional space to maneuver the Legpress athletic.







Adapted lepress atheltic beside the module Back Flexion/Extension and the Back rotation as a full set.





Legpress athletic in Rehabilitation clinic in Dubai with national soccer player







Legpress athletic in Berlin Olympic Trainings center with athlets testing



## Legpress athletic in Rehabilitation clinic with national alpin ski driver





## Legpress athletic in Olympic – Center in Berlin with athlete hurdle runner

